

Serie E Drying oven / Hot air sterilizer E 28

Mit dem Brutschrank, Serie E 28, erhalten Sie ein zuverlässiges und leistungsfähiges Gerät in solider Ausstattung zu einem sehr günstigen Preis.



Leistungsmerkmale/Ausstattung:

- Temperature range 60 °C (140 °F) up to 230 °C (446 °F)
- Robust, space-saving small unit with hydraulic-mechanical control and adjustable ventilation
- Hydraulic-mechanical thermostat
- Optional independent temperature safety device class 1 (DIN 12880)
- Timer 0 to 120 min
- 2 chrome plated racks included
- Units are stackable
- BINDER test certificate

	E 28
Exterior dimensions	
Width (mm / inch)	580 / 22.8
Height (inclusive feet / castors) (mm / inch)	402 / 15.8
Depth (mm / inch)	425 / 16.7
Plus door handle, I-panel and exhaust duct (mm/inch)	85 / 3.4
Number of doors	1
Interior dimensions	
Width (mm / inch)	400 / 15.8
Height (mm / inch)	280 / 11.0
Depth (mm / inch)	250 / 9.8
Interior volume (l / cu.ft.)	28 / 0.99
Racks, chrome-plated (number standard / max.)	2 / 4
Weight of the unit (empty) (kg / lbs.)	22 / 48.6
Temperature data	
Temperature range, 5 °C/9 °F above ambient up to (°C/°F)	60 - 230
Temperature variation 1) at 150 °C (302 °F) (± °C)	3
Temperature fluctuation (± °C)	1.5
Air change at 70 °C (158 °F) (x/h)	30
Electrical data	
Housing protection acc. to EN 60529	IP 20
Nominal voltage (± 10 %) 50/60 Hz (V)	230
Nominal power (W)	800

1) value without window

All technical specification are specified for units with standard equipment at an ambient temperature of 25 °C (77 °F) and a voltage fluctuation of ±10 %. The temperature data are determined in accordance to factory standard following DIN 12880, part 2 respecting the recommended wall clearances of 10 % of the height, width and depth of the inner chamber. All indications are average values, typical for units produced in series.

We reserve the right to alter technical specifications at all times.

O Option - not available

Technical specifications subject to change